



Doc. No.: QC_4.14.2. Rev. Date: 01.04.2019

SparMED ApS



VAT no.: DK 30 89 85 75
Ryttermarken 2, 3520, Farum, Denmark
Phone: +45 39 40 25 03
e-mail: info@sparmed.dk - www.sparmed.dk

Certificate of Analysis

Product ID: Oosafe ® 35 mm Dish, Non-Treated

Product code: OOPW-TF06

LOT No.: 08226

Expiry date: 2027-02

Storage Conditions: 20-30°C, dry room, no exposure to sunlight

Quality Assurance:

Analysis	Reference No.	Results	Completion date
Proven non-embryotoxic by Mouse Embryo Assay Test. Over 80% embryo development to the expanded blastocyst stage within 96 hours.	SPAR-7943-0322	Passed	27-03-2022
Proved stable human sperm motility: ≥75% sperm motility after 24 hours proven.	S5735-0322SPAR	Passed	25-03-2022
Proven non-toxic by Limulus Amebocyte Lysate (LAL) test. Pass criteria <0.03 EU/device.	E2111-0322SPAR	Passed	24-03-2022
Sterilization by gamma irradiation. Delivered irradiation dose: 8.6 kGy-9.5 kGy. Specified irradiation dose: 8.0 kGy-10.0 kGy.	N/A	Passed	N/A

Quality Control according to ISO 13485:2016

GOosafe with SparMED!

Zahraa Baker

Quality Control Department

SparMED ApS

Date: 03-08-2022

