




Doc. No.: QC\_4.14.2. Rev. Date: 01.04.2019

ApS



sparMED  
GOosafe in IVF Laboratories

SparMED

VAT no.: DK 30 89 85 75  
Ryttermarken 2, 3520, Farum, Denmark  
Phone: +45 39 40 25 03  
e-mail: [info@sparmed.dk](mailto:info@sparmed.dk) - [www.sparmed.dk](http://www.sparmed.dk)

## Certificate of Analysis

Product ID: OOSAFE® Center Well Dish With 2 Compartments

Product code: OOPW-CW03

LOT No.: 08237

Expiry date: 2027-04

Storage Conditions: 20-30°C, dry room, no exposure to sunlight

Quality Assurance:

Analysis	Reference No.	Results	Completion date
Proven non-embryotoxic by Mouse Embryo Assay Test. Over 80% embryo development to the expanded blastocyst stage within 96 hours.	SPAR-8259-0522	Passed	20-05-2022
Proved stable human sperm motility: ≥75% sperm motility after 24 hours proven.	S5824-0522SPAR	Passed	18-05-2022
Proven non-toxic by Limulus Amebocyte Lysate (LAL) test. Pass criteria <0.03 EU/device.	E2253-0522SPAR	Passed	17-05-2022
Sterilization by gamma irradiation. Delivered irradiation dose: 8.6 kGy-9.5 kGy. Specified irradiation dose: 8.0 kGy-10.0 kGy.	N/A	Passed	N/A

Quality Control according to ISO 13485:2016

**GOosafe with SparMED!**

Zahraa Baker  
Quality Control Department  
SparMED ApS  
Date: 08-07-2022

