



Doc. No.: QC\_4.14.2. Rev. No. 02 Date: 22.06.2022



sparMED

GOosafe in IVF Laboratories

SparMED ApS

VAT no.: DK 30 89 85 75

Ryttermarken 2, 3520, Farum, Denmark

Phone: +45 39 40 25 03

e-mail: [info@sparmed.dk](mailto:info@sparmed.dk) - [www.sparmed.dk](http://www.sparmed.dk)

## Certificate of Analysis

**Product ID:** Oosafe® Plasticware: Center well dish with 2 compartments

**Product code:** OOPW-CW03

**LOT No.:**08262

**Expiry date:** 2027-07

**Storage Conditions:** 15-30°C, dry room, no exposure to sunlight

**Quality Assurance:**

Analysis	Reference No.	Results	Completion date
Proven non-embryotoxic by Mouse Embryo Assay Test. Over 80% embryo development to the expanded blastocyst stage within 96 hours.	SPAR-9289-0822	Passed	14-08-2022
Proven stable human sperm motility: ≥75% sperm motility after 24 hours proven.	S6012-0822SPAR	Passed	11-08-2022
Proven non-toxic by Limulus Amebocyte Lysate (LAL) test. Pass criteria <0.03 EU/device.	E2503-0822SPAR	Passed	11-08-2022
Sterilized by E-Beam sterilization. Specified sterilization dose: 10.0 kGy-12.0 kGy.	N/A	Passed	N/A

Quality Control according to ISO 13485:2016

**GOosafe with SparMED!**

Andrius Daugintis,

Quality Control

SparMED ApS

Date: 02/09/2022

