

Certificate of Analysis

Product ID: Oosafe® Plasticware: Center well dish with 2 compartments, label area grip

Product code: OOPW-CW05

LOT No.:08263

Expiry date: 2027-07

Storage Conditions: 15-30°C, dry room, no exposure to sunlight

Quality Assurance:

| Analysis | Reference No. | Results | Completion date |
|--|----------------|---------|-----------------|
| Proven non-embryotoxic by Mouse Embryo Assay Test. Over 80% embryo development to the expanded blastocyst stage within 96 hours. | SPAR-9289-0822 | Passed | 14-08-2022 |
| Proven stable human sperm motility: ≥75% sperm motility after 24 hours proven. | S6012-0822SPAR | Passed | 11-08-2022 |
| Proven non-toxic by Limulus Amebocyte Lysate (LAL) test. Pass criteria <0.03 EU/device. | E2503-0822SPAR | Passed | 11-08-2022 |
| Sterilized by E-Beam sterilization. Specified sterilization dose: 10.0 kGy-12.0 kGy. | N/A | Passed | N/A |

Quality Control according to ISO 13485:2016

GOosafe with SparMED!

Andrius Daugintis,

Quality Control

SparMED ApS

Date: 02/09/2022

