

Certificate of Analysis

Product name: 60mm Dish, Non- Treated

Product type: Oosafe® Plasticware IVD

Product code: OOPW-ST06

Batch/LOT No.: 09005

Product expiry date: 2029-06

| TESTING METHOD | ACCEPTANCE CRITERIA | RESULTS | CONCLUSION | DOCUMENT NO. |
|-----------------------------------|--|---------|-----------------------------|----------------|
| Mouse Embryo Assay Test (MEA) | Over 80% embryo development to the expanded blastocyst stage within 96 hours | Passed | Non-embryotoxic | SPAR-4670-0724 |
| Human Sperm Survival Assay (HSSA) | ≥75% sperm motility after 24 hours proven | Passed | Stable human sperm motility | S8118-0724SPAR |
| Limulus Amoebocyte Lysate (LAL) | <0.03 EU/device | Passed | Do not contain endotoxin | E4124-0724SPAR |

Quality Management system according to ISO 13485:2016

Sterilized by irradiation according to ISO 11137-1

Sterilization method: E-Beam

Sterility assurance level: 10⁻⁶


GOosafe with SparMED!

Manuela Jagieniak

Quality Control Department

SparMED ApS

Issue date: 31.07.2024



Sparmed ApS

Ryttermarken 2, 3520 Farum, DK +45 39 40 25 03 info@sparmed.dk www.sparmed.dk
Corp. VAT no.: DK 30 89 85 75