

Certificate of Analysis

Product name: Snap-Cup Tube, 14 ml
Product type: Oosafe® Plasticware IVD
Product code: OOPW-OT10
Batch/LOT No.: 09035
Product expiry date: 2029-10

TESTING METHOD	ACCEPTANCE CRITERIA	RESULTS	CONCLUSION	DOCUMENT NO.
Mouse Embryo Assay Test (MEA)	Over 80% embryo development to the expanded blastocyst stage within 96 hours	Passed	Non-embryotoxic	SPAR-5130-1024
Human Sperm Survival Assay (HSSA)	≥75% sperm motility after 24 hours proven	Passed	Stable human sperm motility	S8289-1024SPAR
Limulus Amoebocyte Lysate (LAL)	<0.03 EU/device	Passed	Do not contain endotoxin	E4290-1024SPAR

Quality Management system according to ISO 13485:2016

Sterilized by irradiation according to ISO 11137-1

Sterilization method: E-Beam

Sterility assurance level: 10⁻⁶


GOosafe with SparMED!

Manuela Jagieniak

Quality Control Department

SparMED ApS

Issue date: 29.10.2024



Sparmed ApS

Ryttermarken 2, 3520 Farum, DK ● +45 39 40 25 03 ● info@sparmed.dk ● www.sparmed.dk
 Corp. VAT no.: DK 30 89 85 75