

# Certificate of Analysis

**Product name:** Snap-Cup Tube, 14 ml  
**Product type:** Oosafe® Plasticware  
**Product code:** OOPW-OT10  
**Batch/LOT No.:** 09201  
**Product expiry date:** 2030-11

TESTING METHOD	ACCEPTANCE CRITERIA	RESULTS	CONCLUSION	TEST REPORT NO.
Mouse Embryo Assay (MEA)	Over 80% embryo development to the expanded blastocyst stage within 96 hours	Passed	Non-embryotoxic	MEA.025.2560.2025
Human Sperm Assay (HSA)	≥75% sperm motility after 24 hours proven	Passed	Stable human sperm motility	SSMA.025.2561.2025
Limulus Amoebocyte Lysate (LAL)	≤ 20 EU/device	Passed	Do not contain endotoxin	LAL.025.2562.2025

Quality Management system according to ISO 13485:2016

**Sterilized by irradiation according to ISO 11137-1**  
**Sterilization method:** E-Beam  
**Sterility assurance level:** 10<sup>-6</sup>

**GOosafe with SparMED!**  
 Manuela Jagieniak  
 Quality Control Department  
 SparMED ApS  
 Issue date: 30.03.2026

