

Gynetics Medical Products CERTIFICATE OF QUALITY

ITEM name and code: Embryo Transfer Catheter The Tulip Memo #4000 Memo A

LOT NUMBER: 20170607

STERILIZATION DATE: 10-06-2017

EXPIRY DATE: 07-06-2021

STERILIZATION TECHNIQUE: Gamma Irradiation

STERILIZATION DOSE: 25 - 39 KGray

QUALITY CONTROL TESTING PROCEDURES

TEST	RESULT	CONCLUSION
Endotoxin Assay (LAL)	<0.5 EU/device	Pass

TEST	SPECIFICATIONS	CONCLUSION
Mouse Embryo Assay	> 80%	Pass
MEA Control	> 80%	Pass



Greetje Hermans
Gynetics Medical Products N.V.

07-07-2017
Date