

Gynetics Medical Products CERTIFICATE OF QUALITY

ITEM name and code: #4164 Probet
Endometrial Curette

LOT NUMBER: 20170607

EXPIRY DATE: 07-06-2021

STERILIZATION DATE: 10-06-2017

STERILIZATION TECHNIQUE: EO

**Results of sterility test with bio-
indicators** No growth of organisms

Gynetics guarantees the following results due to a validated EO sterilization process.

QUALITY CONTROL TESTING PROCEDURES

| TEST | SPECIFICATIONS | Result |
|--------------------------------------|----------------|--------|
| Ethylene oxide residue testing | <250 ppm | Pass |
| Ethylene chlorhydrin residue testing | <250 ppm | Pass |
| Suction function testing | <-110 mbars | Pass |



Greetje Hermans
Gynetics Medical Products N.V.

07-07-2017
Date