

## Gynetics Medical Products CERTIFICATE OF QUALITY

**ITEM name and code:** #4164 Probet  
Endometrial Curette

**LOT NUMBER:** 20170920

**EXPIRY DATE:** 20-09-2021

**STERILIZATION DATE:** 23-09-2017

**STERILIZATION TECHNIQUE:** EO

**Results of sterility test with bio-  
indicators** No growth of organisms

Gynetics guarantees the following results due to a validated EO sterilization process.

### QUALITY CONTROL TESTING PROCEDURES

TEST	SPECIFICATIONS	Result
Ethylene oxide residue testing	<250 ppm	Pass
Ethylene chlorhydrin residue testing	<250 ppm	Pass
Suction function testing	<-110 mbars	Pass

  
Greetje Hermans  
Gynetics Medical Products N.V.

17-11-2017

Date