

Gynetics Medical Products CERTIFICATE OF QUALITY

ITEM name and code: Embryo Transfer Catheter #4000 The Tulip Memo A

LOT NUMBER: 20171205

EXPIRY DATE: 5-12-2021

STERILIZATION DATE: 09-12-2017

STERILIZATION TECHNIQUE: Gamma Irradiation

STERILIZATION DOSE 25 - 39 KGray

QUALITY CONTROL TESTING PROCEDURES

TEST	RESULT	CONCLUSION
Endotoxin Assay (LAL)	<0.5 EU/device	Pass

TEST	SPECIFICATIONS	CONCLUSION
Mouse Embryo Assay	> 80%	Pass
MEA Control	> 80%	Pass


Greetje Hermans
Gynetics Medical Products N.V.

16-5-2018
Date