

CERTIFICATE OF QUALITY

ITEM name and code	#2000 Semtrac Set 4.2 - Soft Embryo Transfer Catheter
LOT NUMBER	20180922
EXPIRY DATE	22 September 2022
STERILIZATION DATE	1 October 2018
STERILIZATION TECHNIQUE	Ethylene Oxide
RESULTS OF BI	No growth of organisms

QUALITY CONTROL TESTING PROCEDURES

TEST	SPECIFICATION	CONCLUSION
Endotoxin Assay (LAL)	<0.5 EU/device	Pass
TEST	SPECIFICATIONS	CONCLUSION
Mouse Embryo Assay	> 80%	Pass
MEA Control	> 80%	Pass



Greetje Hermans
Gynétics Medical Products N.V.

Signature date:
8 November 2018