


CERTIFICATE OF QUALITY

| | |
|--------------------------------|--|
| ITEM name and code | #4225 Smooze Long - Intra-Uterine Insemination Catheter/Cannula Soft |
| LOT NUMBER | 20190213 |
| EXPIRY DATE | 13 February 2023 |
| STERILIZATION DATE | 11 March 2019 |
| STERILIZATION TECHNIQUE | Ethylene Oxide |
| RESULTS OF BI | No growth of organisms |

QUALITY CONTROL TESTING PROCEDURES

| TEST | SPECIFICATION | CONCLUSION |
|-----------------------|----------------------|-------------------|
| Endotoxin Assay (LAL) | <0.5 EU/device | Pass |

| TEST | SPECIFICATIONS | CONCLUSION |
|--------------------|-----------------------|-------------------|
| Mouse Embryo Assay | > 80% | Pass |
| MEA Control | > 80% | Pass |



Evelien Sommen
Gynétics Medical Products N.V.

Signature date:
3 April 2019