

CERTIFICATE OF QUALITY

| | |
|--------------------------------|--|
| ITEM name and code | #4164 Probet - Endometrial Suction Curette |
| LOT NUMBER | 20210825 |
| EXPIRY DATE | 25 August 2025 |
| STERILIZATION DATE | 31 August 2021 |
| STERILIZATION TECHNIQUE | Ethylene Oxide |
| | Biological Indicator testing: no growth of organisms |

QUALITY CONTROL TESTING PROCEDURES

| TEST | SPECIFICATION | CONCLUSION |
|-----------------------|----------------|------------|
| Endotoxin Assay (LAL) | <1.0 EU/device | Pass * |

* LAL results are conform the internal quality control procedure in which testing is conducted based on a representative sampling plan



Merve Toprak
 Gynétics Medical Products N.V.

Signature date:
 28 September 2021