



## CERTIFICATE OF ANALYSIS

**PRODUCT DESCRIPTION:** EZ-Tip® (Vial) Denudation and Transfer Pipette

**LOT NUMBER(S):** 203355, 203356, 203357, 203358, 203359, 203360, 203361, 203362, 203363, 203364.

**EXPIRY DATE:** 2021-02-28


**INDICATIONS FOR USE:** The EZ-Tip pipettes are for denudation, i.e. removing the cumulus from an oocyte prior to Intracytoplasmic sperm injection (ICSI) and in vitro fertilization (IVF), and for handling gametes, embryos and biopsied cells (polar bodies, blastomeres and trophoctoderm) during assisted reproductive techniques (ART).

**QUALITY ASSURANCE:** Sterilised by gamma irradiation using dose between 25 and 35kGy to a SAL of  $10^{-6}$

Assay Test. Pass criteria:  
1-cell MEA  $\geq$  80% hatched blastocysts at 120 hrs: **95% PASS**

Limulus Amoebocyte Lysate Test (LAL). Pass criteria  
<20Eu/device. **PASS**

**STORAGE CONDITIONS:** Temperature 10-30°C  
Humidity 20-80% RH  
Keep away from direct sunlight (UV)

Quality Assurance: 

Date: 01 May 2018

**R<sub>X</sub>**Only *Caution: US Federal law restricts this device to sale by or on the order of a physician.*

Research Instruments Ltd, Bickland Industrial Park, Falmouth, Cornwall, TR11 4TA, UK

Tel : +44 1326 372 753 Fax: +44 1326 378 783

[customer.service@origio.com](mailto:customer.service@origio.com)

[www.origio.com](http://www.origio.com)



CooperSurgical Fertility Companies