



SUCCESS THROUGH INNOVATION

CERTIFICATE OF ANALYSIS

PRODUCT DESCRIPTION:	EZ-Tip (Singles) Denudation and Transfer Pipette 170µm
REF:	7-72-4170
INDICATIONS FOR USE:	The EZ-Tip pipettes are for denudation, i.e. removing the cumulus from an oocyte prior to the assisted reproduction techniques of intracytoplasmic sperm injection (ICSI) and in vitro fertilization (IVF), and effecting dish to dish transfer of cells, embryos and oocytes.
LOT NUMBER:	79996
EXPIRY DATE:	2019-05
QUALITY ASSURANCE:	<p>Sterilised by gamma irradiation using dose between 25 and 35kGy.</p> <p>Pipette tips non-embryotoxic proven by Mouse Embryo Assay Test (RI test protocol ref RIMEA1). Pass criteria 2-cells at day 2 \geq 80%, expanded blastocysts at day 6 \geq 80%. PASS</p> <p>Pipette tips non-pyrogenic proven by Limulus Amoebocyte Lysate Test. Pass criteria $<$20Eu/device. PASS.</p>
STORAGE CONDITIONS:	Temperature 10-30°C Humidity 30-60% RH Keep away from direct sunlight (UV)

Quality Assurance:

Date: 13th June 2014

R_xOnly *Caution: US Federal law restricts this device to sale by or on the order of a physician.*

Research Instruments Ltd, Bickland Industrial Park, Falmouth, Cornwall, TR11 4TA, UK

Tel : +44 1326 372 753

Fax: +44 1326 378 783

sales@research-instruments.comwww.research-instruments.com