



## Certificate of Analysis

### Product: SpermCryo™ All-round

Batch number: CRY103  
Expiration Date: 31-10-2022  
(DD/MM/YYYY)

| Analysis                              | Limits       | Results     |
|---------------------------------------|--------------|-------------|
| pH:                                   | 7,20 – 7,60  | 7,53        |
| Sterility Test Bacteria:              | NEGATIVE     | PASSED      |
| Sterility Test Fungi:                 | NEGATIVE     | PASSED      |
| Endotoxin:                            | < 0,25 EU/ml | < 0,1 EU/ml |
| Sperm Survival:<br>(4 hours exposure) | ≥ 80%        | 95%         |

---

Approved on: 29-04-2021  
(DD/MM/YYYY)  
By: M. Berbers

Signature:

A handwritten signature in blue ink, appearing to read 'M. Berbers', is written over a faint, rectangular stamp.

Quality Control Department, Gynotec B.V.