



Certificate of Analysis

Product: SpermCryo[®] All-round

Batch number: CRY2402
Expiration Date: 31-08-2025
(DD/MM/YYYY)

| Analysis | Limits | Results |
|--|---------------------|---------|
| pH: | 7,20 – 7,60 | 7,49 |
| Sterility Test (According to current Ph. Eur. 2.6.1 / USP <71>): | No Growth | PASSED |
| Endotoxin: | < 0,25 EU/ml | < 0,10 |
| HSA concentration: | 3,20 – 4,80 g/liter | 4,00 |
| HSSA (<i>motility compared with control after 4h exposure to test medium</i>): | ≥ 80% | 100% |

Approved on: 13-08-2024
(DD/MM/YYYY)

By: M. Berbers

Signature:

A handwritten signature in blue ink, appearing to read 'M. Berbers', is written over a faint, circular stamp or watermark.