



## Certificate of Analysis

### Product: SpermFilter® Gradient 45%

Batch number: FLT45-223  
Expiration Date: 31-01-2024  
(DD/MM/YYYY)

Analysis	Limits	Results
pH:	7,20 – 7,60	7,41
Osmolality:	310 – 340 mOsm/kg	325
Sterility Test Bacteria:	NEGATIVE	PASSED
Sterility Test Fungi:	NEGATIVE	PASSED
<i>(Sterile by filtration over 0,2µm filter)</i>		
Endotoxin:	< 0,5 EU/ml	< 0,100 EU/ml
Sperm Survival:	≥ 80%	98%
<i>(4 hours exposure)</i>		
Sperm Survival:	≥ 75%	95%
<i>(24 hours exposure)</i>		

---

Approved on: 28-07-2022  
(DD/MM/YYYY)  
By: M. Berbers

Signature:

Quality Control Department, Gynotec B.V.