



## Certificate of Analysis

### Product: SpermFilter® Gradient 80%

Batch number: FLT80-010  
Expiration Date: 30-09-2021  
(DD/MM/YYYY)

Analysis	Limits	Results
pH:	7,20 – 7,60	7,34
Osmolality:	310 – 330 mOsm/kg	324
Sterility Test Bacteria:	NEGATIVE	PASSED
Sterility Test Fungi:	NEGATIVE	PASSED
Endotoxin:	< 0,5 EU/ml	< 0,100 EU/ml
Density:	1,0970 – 1,1070 g/ml	1,1025
Sperm Survival:	≥ 80%	98%

---

Approved on: 31-03-2020  
(DD/MM/YYYY)

By: J. Berbers

Signature:

Quality Control Department, Gynotec B.V.