



## Certificate of Analysis

### Product: SpermFilter® Gradient 80%

Batch number: FLT80-021  
Expiration Date: 31-01-2022  
(DD/MM/YYYY)

| Analysis                 | Limits               | Results       |
|--------------------------|----------------------|---------------|
| pH:                      | 7,20 – 7,60          | 7,34          |
| Osmolality:              | 310 – 330 mOsm/kg    | 317           |
| Sterility Test Bacteria: | NEGATIVE             | PASSED        |
| Sterility Test Fungi:    | NEGATIVE             | PASSED        |
| Endotoxin:               | < 0,5 EU/ml          | < 0,100 EU/ml |
| Density:                 | 1,0970 – 1,1070 g/ml | 1,0986        |
| Sperm Survival:          | ≥ 80%                | 99%           |

---

Approved on: 31-07-2020  
(DD/MM/YYYY)

By: J. Berbers

Signature:

A handwritten signature in blue ink, appearing to read 'J. Berbers', is written over a light blue horizontal line.

Quality Control Department, Gynotec B.V.