



## Certificate of Analysis

### Product: SpermFilter® Gradient 80%

Batch number: FLT80-111  
Expiration Date: 30-09-2022  
(DD/MM/YYYY)

Analysis	Limits	Results
pH:	7,20 – 7,60	7,33
Osmolality:	310 – 330 mOsm/kg	328
Sterility Test Bacteria:	NEGATIVE	PASSED
Sterility Test Fungi:	NEGATIVE	PASSED
Endotoxin:	< 0,5 EU/ml	< 0,100 EU/ml
Density:	1,0970 – 1,1070 g/ml	1,1023
Sperm Survival:	≥ 80%	99%

---

Approved on: 22-03-2021  
(DD/MM/YYYY)

By: M. Berbers

Signature:

Quality Control Department, Gynotec B.V.