



Certificate of Analysis

Product: SpermFilter[®] Gradient 80%

Batch number: FLT80-139
Expiration Date: 30-06-2023
(DD/MM/YYYY)

Analysis	Limits	Results
pH:	7,20 – 7,60	7,37
Osmolality:	310 – 340 mOsm/kg	323
Sterility Test Bacteria:	NEGATIVE	PASSED
Sterility Test Fungi:	NEGATIVE	PASSED
(Sterile by filtration over 0,2µm filter)		
Endotoxin:	< 0,5 EU/ml	< 0,100 EU/ml
Density:	1,0970 – 1,1070 g/ml	1,1032
Viscosity in cP (25°)	< 1,65 cP	1,50
Sperm Survival:	≥ 80%	97%
(4 hours exposure)		

Approved on: 20-12-2021
(DD/MM/YYYY)

By: M. Berbers

Signature:

Quality Control Department, Gynotec B.V.