



## Certificate of Analysis

### Product: SpermFilter<sup>®</sup> Gradient 80%

Batch number: FLT80-205  
Expiration Date: 31-08-2023  
(DD/MM/YYYY)

Analysis	Limits	Results
pH:	7,20 – 7,60	7,34
Osmolality:	310 – 340 mOsm/kg	326
Sterility Test Bacteria:	NEGATIVE	PASSED
Sterility Test Fungi:	NEGATIVE	PASSED
<i>(Sterile by filtration over 0,2µm filter)</i>		
Endotoxin:	< 0,5 EU/ml	< 0,100 EU/ml
Density:	1,0970 – 1,1070 g/ml	1,1001
Viscosity in cP (25°)	< 1,65 cP	1,53
Sperm Survival:	≥ 80%	100%
<i>(4 hours exposure)</i>		

---

Approved on: 14-03-2022  
(DD/MM/YYYY)

By: M. Berbers

Signature:

Quality Control Department, Gynotec B.V.