



Certificate of Analysis

Product: SpermFilter® Gradient 80%

Batch number:

FLT80-824

Expiration Date:
(DD/MM/YYYY)

30-11-2019

Analysis	Limits	Results
pH:	7,20 – 7,60	7,39
Osmolality:	310 – 330 mOsm/kg	322
Sterility Test Bacteria:	NEGATIVE	PASSED
Sterility Test Fungi:	NEGATIVE	PASSED
Endotoxin:	< 0,5 EU/ml	< 0,1 EU/ml
Density:	1,0970 – 1,1070 g/ml	1,0996
Sperm Survival:	≥ 80%	99%

Approved on:
(DD/MM/YYYY)

07-06-2018

By:

J. Berbers

Signature:

Quality Control Department, Gynotec B. V.