

PRODUCT: Vitrification Media

LOT#: G0103

EXPIRY DATE: 2017.04.26

Year / Month / Date

REF: 91171

CODE: VT801



(01) 1 4582231 46076 9 (17) 170426 (10) G0103

APPROVED BY:

Eri Nishimura
Quality Assurance

2017.01.26

CERTIFICATE OF ANALYSIS

Indication for Use: Vitrification Media is used for
vitrification and storage of oocytes and/or embryos.

Storage: Store between 2-8 °C.

Caution: Federal law restricts this device to sale by or on the order of a physician or practitioner
trained in its use.

We certify that products manufactured by Kitazato BioPharma Co., Ltd. have been approved in accordance with current Good Manufacturing Practice (cGMP). The quality of products was tested and evaluated by quality control procedures.

| ASSAY | SPECIFICATION | RESULT | | | |
|----------------------------|------------------------------------|--------|--------|--------|--------|
| | | BS | ES | VS1 | VS2 |
| Color | Clear, Particle Free/ Rose Pink | Pass | Pass | Pass | Pass |
| pH | 7.20-7.60 | 7.41 | 7.40 | 7.48 | 7.49 |
| Osmolarity | BS:240-300 | 262 | 2706 | 5583 | 5700 |
| | ES:2300-2800 VS:4900-6000 | | | | |
| Endotoxin ¹ | < 0.25 EU/mL | <0.006 | <0.006 | <0.006 | <0.006 |
| Sterility ² | No Microbial Growth | Pass | Pass | Pass | Pass |
| MEA(One Cell) ³ | | | | | |
| -Control | ≥ 80% | | 100% | | |
| -This Lot | ≥ 80% | | 100% | | |

¹ Bacteria Endotoxins Test, Photometric Techniques in accordance with the USP, EP and JP

² Sterility Tests, Direct Inoculation in accordance with the USP, EP and JP

³ One cell mouse embryo testing is performed using fresh mouse embryos. Test results indicate the percentage of mouse embryos developing to fully expanded blastocysts after 96 hours in culture medium.