

**Certificate of Analysis
GM501 Mineral Oil**

Batchnumber OIL-1737
 Production date 11 December 2017
 Expiry date 30 June 2019
 Total volume (ml) 52500

	Results	Specifications
Appearance	PASSED	PASSED
Sterility test Sterilised by sterile filtration over 0,2µm filter	PASSED	PASSED
Acidity	Conform Ph Eur	Conform Ph Eur
Solid paraffin	Conform Ph Eur	Conform Ph Eur
Viscosity	20,2	< 30 cP at 30°C
Specific gravity	0,8514	< 0,86 g/ml
Endotoxin	< 0,05	< 0,25 EU/ml
Peroxide value (POV)	< 0,1	< 0,1 mEq/kg
Mouse Embryo Assay (1-cell to expanded blastocyst after 120 hours)	80,95	≥ 80 %

Release date 27 December 2017

Released by Vanessa Bierlich

Signature