

**Certificate of Analysis
GM501 Mineral Oil**

Batchnumber	OIL-1738
Production date	14 December 2017
Expiry date	30 June 2019
Total volume (ml)	150000

	Results	Specifications
Appearance	PASSED	PASSED
Sterility test <small>Sterilised by sterile filtration over 0,2µm filter</small>	PASSED	PASSED
Acidity	Conform Ph Eur	Conform Ph Eur
Solid paraffin	Conform Ph Eur	Conform Ph Eur
Viscosity	21,2	< 30 cP at 30°C
Specific gravity	0,851	< 0,86 g/ml
Endotoxin	< 0,1	< 0,25 EU/ml
Peroxide value (POV)	< 0,1	< 0,1 mEq/kg
Mouse Embryo Assay <small>(1-cell to expanded blastocyst after 120 hours)</small>	85,71	≥ 80 %

Release date 03 January 2018

Released by Vanessa Bierlich

Signature