



## Certificate of Analysis

### Product: SpermWash®

Batch number: SW-2462  
Expiration Date: 31-03-2026  
(DD/MM/YYYY)

| Analysis                                                                     | Limits              | Results |
|------------------------------------------------------------------------------|---------------------|---------|
| pH:                                                                          | 7,30 – 7,60         | 7,44    |
| Osmolality:                                                                  | 270 – 290 mOsm/kg   | 280     |
| Sterility Test (According to current Ph. Eur. 2.6.1 / USP <71>):             | No Growth           | PASSED  |
| Endotoxin:                                                                   | < 0,25 EU/ml        | < 0,05  |
| HSA concentration:                                                           | 3,20 – 4,80 g/liter | 3,80    |
| Gentamicin concentration:                                                    | 5 – 15 mg/liter     | 10,05   |
| MEA:<br>(% blastocysts after 96 hours; 60 min. exposure time to test medium) | ≥ 80%               | 81%     |
| HSSA (motility compared with control after 24h exposure to test medium):     | ≥ 80%               | 100%    |

---

Approved on: 11-10-2024  
(DD/MM/YYYY)

By: M. Berbers

Signature:

A handwritten signature in blue ink, appearing to read 'M. Berbers', is written over a faint, circular stamp or watermark.