



Certificate of Analysis

Product: SpermWash®

Batch number: SW-2475
Expiration Date: 30-06-2026
(DD/MM/YYYY)

| Analysis | Limits | Results |
|--|---------------------|---------|
| pH: | 7,30 – 7,60 | 7,44 |
| Osmolality: | 270 – 290 mOsm/kg | 278 |
| Sterility Test (According to current Ph. Eur. 2.6.1 / USP <71>): | No Growth | PASSED |
| Endotoxin: | < 0,25 EU/ml | < 0,05 |
| HSA concentration: | 3,20 – 4,80 g/liter | 3,90 |
| Gentamicin concentration: | 5 – 15 mg/liter | 12,80 |
| MEA: (% blastocysts after 96 hours; 60 min. exposure time to test medium) | ≥ 80% | 100% |
| HSSA (motility compared with control after 24h exposure to test medium): | ≥ 80% | 98% |

Approved on: 18-12-2024
(DD/MM/YYYY)

By: M. Berbers

Signature:

A handwritten signature in blue ink, appearing to read 'M. Berbers', is written over a faint, circular stamp or watermark.