



## Certificate of Analysis

### Product: SpermWash®

Batch number: SW-2534  
Expiration Date: 30-11-2026  
(DD/MM/YYYY)

| Analysis                                                                 | Limits              | Results |
|--------------------------------------------------------------------------|---------------------|---------|
| pH:                                                                      | 7,30 – 7,60         | 7,41    |
| Osmolality:                                                              | 270 – 290 mOsm/kg   | 281     |
| Sterility Test (According to current Ph. Eur. 2.6.1 / USP <71>):         | No Growth           | PASSED  |
| Endotoxin:                                                               | < 0,25 EU/ml        | < 0,05  |
| HSA concentration:                                                       | 3,20 – 4,80 g/liter | 3,90    |
| Gentamicin concentration:                                                | 5 – 15 mg/liter     | 9,00    |
| MEA:                                                                     | ≥ 80%               | 100%    |
| (% blastocysts after 96 hours; 60 min. exposure time to test medium)     |                     |         |
| HSSA (motility compared with control after 24h exposure to test medium): | ≥ 80%               | 91%     |

---

Approved on: 10-06-2025  
(DD/MM/YYYY)  
By: M. Berbers

Signature: