



## Certificate of Analysis

### Product: SpermWash®

Batch number: SW-2624  
Expiration Date: 31-10-2027  
(DD/MM/YYYY)

| Analysis  | Limits              | Results |
|---|---------------------|---------|
| pH:   | 7,30 – 7,60         | 7,43    |
| Osmolality:   | 270 – 290 mOsm/kg   | 281     |
| Sterility Test (According to current Ph. Eur. 2.6.1 / USP <71>):                  | No Growth           | PASSED  |
| Endotoxin:  | < 0,25 EU/ml        | < 0,05  |
| HSA concentration:  | 3,20 – 4,80 g/liter | 3,90    |
| Gentamicin concentration:   | 5 – 15 mg/liter     | 11,00   |
| MEA:<br>(% blastocysts after 96 hours; 60 min. exposure time to test medium)      | ≥ 80%               | 100%    |
| HSSA ( <i>motility compared with control after 24h exposure to test medium</i> ): | ≥ 80%               | 100%    |

---

Approved on: 06-05-2026  
(DD/MM/YYYY)  
By: M. Berbers

Signature: