



## Certificate of Analysis

### Product: SpermWash<sup>®</sup>

Batch number: SW-334  
Expiration Date: 31-12-2024  
(DD/MM/YYYY)

| Analysis                                                                        | Limits              | Results |
|---------------------------------------------------------------------------------|---------------------|---------|
| pH:                                                                             | 7,30 – 7,60         | 7,49    |
| Osmolality:                                                                     | 270 – 290 mOsm/kg   | 284     |
| Sterility Test (According to<br>current Ph. Eur. 2.6.1 / USP <71>):             | No Growth           | PASSED  |
| Endotoxin:                                                                      | < 0,25 EU/ml        | < 0,05  |
| HSA concentration:                                                              | 3,20 – 4,80 g/liter | 4,20    |
| MEA:<br>(% blastocysts after 96 hours; 60<br>min. exposure time to test medium) | ≥ 80%               | 100%    |

---

Approved on: 28-06-2023  
(DD/MM/YYYY)  
By: M. Berbers

Signature: