

## Certificate of Analysis

**Product name: OOSAFE® DISINFECTION WIPES**

**Product code: OODW-70**

**Batch/LOT No.: 01803**

**Production date: 2020-05**

**Expiry date: 2022-09**

Analysis / Method	Specification / Reference No.	Result	Completion date
Appearance / Visual	Clear liquid inside wipes	Clear liquid inside wipes	NA
Odor / Olfactive	Odorless	Odorless	NA
pH-Rate (2% active substance in water) / pH meter (DIN EN 1262)	7,0-9,0	8,3	NA
Density, 20 °C, gr/cm <sup>3</sup> / Density measuring device	1,001-1,010	1,004	NA
Quaternary Ammonium Compounds / QUAT test stripe	1,9-2,2 %	2 %	NA
Water		the rest	NA
Proven non-embryotoxic by Mouse Embryo Assay Test. Over 80% embryo development to the expanded blastocyst stage within 96 hours.	SPAR-4009-0520	Passed	19-05-2020
Proved stable human sperm motility: ≥75% sperm motility after 24 hours proven.	S4456-0520SPAR	Passed	20-05-2020

Quality Control according to ISO 13458:2016

**GOosafe with SparMED!**

Inga Sileikaite,  
Quality Control Department  
SparMED ApS  
**Date: 21-May-2020**

