

Certificate of Analysis

Product name: OOSAFE® DISINFECTION WIPES**Product code: OODW-70****Batch/LOT No.: 03808****Production date: 2020-09****Expiry date: 2023-09**

| Analysis / Method | Specification / Reference No. | Result | Completion date |
|--|-------------------------------|---------------------------|-----------------|
| Appearance / Visual | Clear liquid inside wipes | Clear liquid inside wipes | NA |
| Odor / Olfactive | Odorless | Odorless | NA |
| pH-Rate (2% active substance in water) / pH meter (DIN EN 1262) | 7,0 – 9,0 | 8,2 | NA |
| Density, 20 °C, gr/cm ³ / Density measuring device | 1,001 – 1,010 | 1,004 | NA |
| Quaternary Ammonium Compounds / QUAT test stripe | 1,9 – 2,2 % | 2 % | NA |
| Water | | the rest | NA |
| Proven non-embryotoxic by Mouse Embryo Assay Test. Over 80% embryo development to the expanded blastocyst stage within 96 hours. | SPAR-4567-0920 | Passed | 03-10-2020 |
| Proved stable human sperm motility: ≥75% sperm motility after 24 hours proven. | S4633-0920SPAR | Passed | 30-09-2020 |

Quality Control according to ISO 13458:2016

GOosafe with SparMED!Florenca Garcia
Quality Control Department
SparMED ApS
Date: 05-10-2020