

## Certificate of Analysis

**Product name: OOSAFE® DISINFECTANT- WIPES**

**Product code: OOMDW-120**

**Batch/LOT No.: 03809-MD**

**Production date: 2020-09**

**Expiry date: 2023-09**

| Analysis / Method  | Specification / Reference No. | Result                    | Completion date |
|--|-------------------------------|---------------------------|-----------------|
| Appearance / Visual  | Clear liquid inside wipes     | Clear liquid inside wipes | NA              |
| Odor / Olfactive   | Odorless                      | Odorless                  | NA              |
| pH-Rate (2% active substance in water) / pH meter (DIN EN 1262)  | 7,0 – 9,0                     | 8,2                       | NA              |
| Density, 20 °C, gr/cm <sup>3</sup> / Density measuring device  | 1,001 – 1,010                 | 1,004                     | NA              |
| Quaternary Ammonium Compounds / QUAT test stripe   | 1,9 – 2,2 %                   | 2 %                       | NA              |
| Water  |                               | the rest                  | NA              |
| Proven non-embryotoxic by Mouse Embryo Assay Test. Over 80% embryo development to the expanded blastocyst stage within 96 hours. | SPAR-4567-0920                | Passed                    | 03-10-2020      |
| Proved stable human sperm motility: ≥75% sperm motility after 24 hours proven.   | S4633-0920SPAR                | Passed                    | 30-09-2020      |

Quality Control according to ISO 13458:2016

### GOosafe with SparMED!

Florencia Garcia  
Quality Control Department  
SparMED ApS  
Date: 05-10-2020


