

Certificate of Analysis

Product name: Disinfection wipes**Product code: OODW-70****Batch/LOT No.: 13603****Production date: 2021-09****Expiry date: 2024-09**

Analysis / Method	Specification / Reference No.	Result	Completion date
Appearance / Visual	Clear liquid inside wipes	Clear liquid inside wipes	15-09-2021
Odor / Olfactive	Odourless	Odourless	15-09-2021
pH-Rate (2% active substance in water) / pH meter (DIN EN 1262)	7.0 – 9.0	8.1	15-09-2021
Density, 20 °C, gr/cm ³ / Density measuring device	1.001 – 1.010	1.008	15-09-2021
Quaternary Ammonium Compounds / QUAT test stripe	1.9 – 2.2 %	2 %	15-09-2021
Water		The rest	
Proven non-embryotoxic by Mouse Embryo Assay Test. Over 80% embryo development to the expanded blastocyst stage within 96 hours.	SPAR-0781-0921	Passed	27-09-2021
Proved stable human sperm motility: ≥75% sperm motility after 24 hours proven.	S5349-0921SPAR	Passed	24-09-2021

Quality Control according to ISO 13458:2016

GOosafe with SparMED!

Andrius Daugintis
Quality Control Department
SparMED ApS
Date: 01-10-2021

