

Certificate of Analysis

Product name: OOMDW-120

Product code: Oosafe® Disinfectant-Wipes, 70 wipes/bucket

Batch/LOT No.: 31802-MD

Production date: 2023-04

Expiry date: 2026-04

Analysis / Method	Specification / Reference No.	Result	Completion date
Appearance / Visual	Clear liquid inside wipes	Clear liquid inside wipes	28/04/2023
Odor / Olfactive	Odourless	Odourless	28/04/2023
pH-Rate (2% active substance in water) / pH meter (DIN EN 1262)	7,0 – 9,0	8,2	28/04/2023
Density, 20 °C, gr/cm ³ / Density measuring device	1.001 – 1.010	1,008	28/04/2023
Quaternary Ammonium Compounds / QUAT test stripe	1,9 – 2,2 %	2%	28/04/2023
Sodium nitrite		NA	
Additives		NA	
Water		The rest	
Proven non-embryotoxic by Mouse Embryo Assay Test. Over 80% embryo development to the expanded blastocyst stage within 96 hours.	SPAR-2027-0523	Passed	05/06/2023
Proved stable human sperm motility: ≥75% sperm motility after 24 hours proven.	S6555-0623SPAR	Passed	07/06/2023

Quality Control according to ISO 13458:2016

GOosafe with SparMED!

Zahraa Baker

Quality Control Department

SparMED ApS

Date: 12/06/2023




	PREPARED BY:	APPROVED BY:
	QUALITY CONTROL MANAGER	MANAGEMENT REPRESENTATIVE
	Signature 	Signature 