

Certificate of Analysis

Product name: OOSAFE® DISINFECTANT WIPES

Product code: OOMDW-120

Batch/LOT No.: 32403-MD

Production date: 05.2023

Expiry date: 05.2026

Analysis / Method	Specification / Reference No.	Result	Completion date
Appearance / Visual	Clear liquid, inside wipes	Clear liquid, inside wipes	04.07.2023
Odor / Olfactive	Odourless	Odourless	04.07.2023
pH-Rate (2% active substance in water) / pH meter (DIN EN 1262)	7,0 - 9,0	8,2	04.07.2023
Density, 20 °C, gr/cm ³ / Density measuring device	1,001 – 1,010	1,008	04.07.2023
Quaternary Ammonium Compounds / QUAT test stripe	1,9 – 2,2 %	2 %	04.07.2023
Water		The rest	N/A
Proven non-embryotoxic by Mouse Embryo Assay Test. Over 80% embryo development to the expanded blastocyst stage within 96 hours.	SPAR- 2262-0723	Passed	14.07.2023
Proved stable human sperm motility: ≥75% sperm motility after 24 hours proven.	S6622-0723SPAR	Passed	12.07.2023

Quality Control according to ISO 13458:2016

GOosafe with SparMED!

Marzena Gabiga

Quality Control Department

SparMED ApS

Date: 18/07/2023



M. Gabiga

 GOosafe in IVF Laboratories	PREPARED BY:	APPROVED BY:
	QUALITY CONTROL MANAGER	MANAGEMENT REPRESENTATIVE
	Signature <i>livei</i>	Signature <i>Bmbe</i>