

# Certificate of Analysis

**Product name: Oosafe® Disinfectant Wipes**

**Product type: Biocide**

**Product code: OODW-70**

**Batch/LOT No.: 40801**

**Production date: 2024-02**

**Product expiry date: 2027-02**

TESTING METHOD	ACCEPTANCE CRITERIA	RESULTS	DOCUMENT NO.
Color / Visual	Clear liquid	Clear liquid	QC_4.1-0324
Odor / Olfactive	odorless	odorless	PRD/F01
ph (2% active substance in water)/ ph meter	7,0-9,0	8,3	PRD/F01
Density, 20 °C, gr/cm <sup>3</sup> / Density measuring device	1,001-1,010	1,008	PRD/F01

TESTING METHOD	ACCEPTANCE CRITERIA	RESULTS	CONCLUSION	DOCUMENT NO.
Mouse Embryo Assay Test (MEA)	Over 80% embryo development to the expanded blastocyst stage within 96 hours	Passed	Non-embryotoxic	SPAR-3823-0324
Human Sperm Survival Assay (HSSA)	≥75% sperm motility after 24 hours proven	Passed	Stable human sperm motility	S7137-0324SPAR

Quality Management system according to ISO 13485:2016


**GOosafe with SparMED!**

Manuela Jagieniak

Quality Control Department

SparMED ApS

**Issue date: 02.04.2024**



**Sparmed ApS**

Ryttermarken 2, 3520 Farum, DK +45 39 40 25 03 info@sparmed.dk www.sparmed.dk

Corp. VAT no.: DK 30 89 85 75