

Certificate of Analysis

Product name: Oosafe® Disinfectant Wipes

Product type: Biocide

Product code: OODW-70

Batch/LOT No.: 40802

Production date: 2024-02

Product expiry date: 2027-02

| TESTING METHOD | ACCEPTANCE CRITERIA | RESULTS | DOCUMENT NO. |
|---|---------------------|--------------|--------------|
| Color / Visual | Clear liquid | Clear liquid | QC_4.1-0324 |
| Odor / Olfactive | odorless | odorless | PRD/F01 |
| ph (2% active substance in water)/ ph meter | 7,0-9,0 | 8,3 | PRD/F01 |
| Density, 20 °C, gr/cm ³ / Density measuring device | 1,001-1,010 | 1,008 | PRD/F01 |

| TESTING METHOD | ACCEPTANCE CRITERIA | RESULTS | CONCLUSION | DOCUMENT NO. |
|-----------------------------------|--|---------|-----------------------------|----------------|
| Mouse Embryo Assay Test (MEA) | Over 80% embryo development to the expanded blastocyst stage within 96 hours | Passed | Non-embryotoxic | SPAR-3822-0324 |
| Human Sperm Survival Assay (HSSA) | ≥75% sperm motility after 24 hours proven | Passed | Stable human sperm motility | S7137-0324SPAR |

Quality Management system according to ISO 13485:2016


GOosafe with SparMED!

Manuela Jagieniak

Quality Control Department

SparMED ApS

Issue date: 02.04.2024



Sparmed ApS

Ryttermarken 2, 3520 Farum, DK ☎ +45 39 40 25 03 ☎ info@sparmed.dk ☎ www.sparmed.dk

Corp. VAT no.: DK 30 89 85 75