

Certificate of Analysis

Product name: Oosafe® Disinfection wipes

Product type: Biocide

Product code: OODW-70

Batch/LOT No.: 42305

Production date: 2024-05

Product expiry date: 2027-05

TESTING METHOD	ACCEPTANCE CRITERIA	RESULTS	DOCUMENT NO.
Color / Visual	Clear liquid inside wipes	Clear liquid inside wipes	QC_4.1-0624
Odor / Olfactive	odorless	odoless	PRD/F01
ph (2% active substance in water)/ ph meter	7,0 – 9,0	8,1	PRD/F01
Density, 20 °C, gr/cm ³ / Density measuring device	11,001 – 1,010	1,008	PRD/F01

TESTING METHOD	ACCEPTANCE CRITERIA	RESULTS	CONCLUSION	DOCUMENT NO.
Mouse Embryo Assay Test (MEA)	Over 80% embryo development to the expanded blastocyst stage within 96 hours	Passed	Non-embryotoxic	SPAR-4500-0624
Human Sperm Survival Assay (HSSA)	≥75% sperm motility after 24 hours proven	Passed	Stable human sperm motility	S8057-0624SPAR

Quality Management system according to ISO 13485:2016


GOosafe with SparMED!

Manuela Jagieniak

Quality Control Department

SparMED ApS

Issue date: 04.07.2024



Sparmed ApS

Page 1 of 1