

# Certificate of Analysis

**Product name:** Oosafe® Disinfectant Wipes

**Product type:** Biocide

**Product code:** OODW-70

**Batch/LOT No.:** 51002

**Production date:** 03.2025

**Product expiry date:** 03.2027

| TESTING METHOD  | ACCEPTANCE CRITERIA       | RESULTS                   | DOCUMENT NO. |
|---|---------------------------|---------------------------|--------------|
| Color / Visual  | Clear liquid inside wipes | Clear liquid inside wipes | QC_4.1-0325  |
| Odor / Olfactive  | odorless                  | odorless                  | PRD/F01      |
| ph (2% active substance in water)/ ph meter                   | 7.0-9.0                   | 7.7                       | PRD/F01      |
| Density, 20 °C, gr/cm <sup>3</sup> / Density measuring device | 1,001-1,010               | 1.009                     | PRD/F01      |

| TESTING METHOD                    | ACCEPTANCE CRITERIA  | RESULTS | CONCLUSION                  | DOCUMENT NO.   |
|-----------------------------------|--|---------|-----------------------------|----------------|
| Mouse Embryo Assay Test (MEA)     | Over 80% embryo development to the expanded blastocyst stage within 96 hours | Passed  | Non-embryotoxic             | SPAR-718-25-03 |
| Human Sperm Survival Assay (HSSA) | ≥75% sperm motility after 24 hours proven                                    | Passed  | Stable human sperm motility | SPAR-719-25-03 |

Quality Management system according to ISO 13485:2016

**GOosafe with SparMED!**

Manuela Jagieniak

Quality Control Department

SparMED ApS

**Issue date:** 19.05.2025



**Sparmed ApS**

Page 1 of 1

Ryttermarken 2, 3520 Farum, DK +45 39 40 25 03 info@sparmed.dk www.sparmed.dk

Corp. VAT no.: DK 30 89 85 75