

# Certificate of Analysis

**Product name: Oosafe® Medical Device Disinfectant**

**Product type: Class IIa**

**Product code: OOMDW-120**

**Batch/LOT No.: 51903-MD**

**Production date: 2025-05**

**Product expiry date: 2028-05**

TESTING METHOD	ACCEPTANCE CRITERIA	RESULTS	DOCUMENT NO.
Color / Visual	Clear liquid inside wipes	Clear liquid inside wipes	QC_4.1-0625
Odor / Olfactive	odorless	odorless	PRD/01
ph (2% active substance in water)/ ph meter	7.0-9.0	7.6	PRD/01
Density, 20 °C, gr/cm <sup>3</sup> / Density measuring device	1.001-1.010	1.007	PRD/01

TESTING METHOD	ACCEPTANCE CRITERIA	RESULTS	CONCLUSION	DOCUMENT NO.
Mouse Embryo Assay Test (MEA)	Over 80% embryo development to the expanded blastocyst stage within 96 hours	Passed	Non-embryotoxic	SPAR-1061-25-06
Human Sperm Survival Assay (HSSA)	≥75% sperm motility after 24 hours proven	Passed	Stable human sperm motility	SPAR-1062-25-06

Quality Management system according to ISO 13485:2016

**GOosafe with SparMED!**

Manuela Jagieniak

Quality Control Department

SparMED ApS

**Issue date: 20.06.2025**



**Sparmed ApS**

Page 1 of 1