

Certificate of Analysis

Product name: Oosafe® Disinfectant Wipes

Product type: MD Class IIa

Product code: OOMDW-120

Batch/LOT No.: 60502-MD

Production date: 2026-01

Product expiry date: 2028-12

TESTING METHOD	ACCEPTANCE CRITERIA	RESULTS	DOCUMENT NO.
Color / Visual	Clear liquid	Clear Liquid	QC_4.1-0326
Odor / Olfactive	Mild	Mild	F-WI-08-01-4
ph (2% active substance in water)/ ph meter	7.0-8.5	7.7	F-WI-08-01-4
Density, 20 °C, gr/cm ³ / Density measuring device	0.96-1.04	1.006	F-WI-08-01-4

TESTING METHOD	ACCEPTANCE CRITERIA	RESULTS	CONCLUSION	DOCUMENT NO.
Mouse Embryo Assay (MEA)	Over 80% embryo development to the expanded blastocyst stage within 96 hours	Passed	Non-embryotoxic	MEA.025.2910.2025
Human Sperm Assay (HSA)	≥75% sperm motility after 24 hours proven	Passed	Stable human sperm motility	SSMA.025.2911.2025

Quality Management system according to ISO 13485:2016

GOosafe with SparMED!

Marzena Gabiga

Quality Control Department

SparMED ApS

Issue date: 03.03.2026



Sparmed ApS

Page 1 of 1