

Certificate of Analysis

Product name: Oosafe® Disinfectant for Co2 Incubators and Laminar Flow Hoods

Product type: Biocide

Product code: OODIH

Batch/LOT No.: 60965

Production date: 2024-05

Product expiry date: 2027-05

TESTING METHOD	ACCEPTANCE CRITERIA	RESULTS	DOCUMENT NO.
Color / Visual	Clear liquide	Clear liquide	QC_4.1-0524
Odor / Olfactive	mild	mild	PRD-F01
ph (2% active substance in water)/ ph meter	7,0 – 8,5	7,7	PRD-F01
Density, 20 °C, gr/cm ³ / Density measuring device	1,01 – 2,8	1,08	PRD-F01

TESTING METHOD	ACCEPTANCE CRITERIA	RESULTS	CONCLUSION	DOCUMENT NO.
Mouse Embryo Assay Test (MEA)	Over 80% embryo development to the expanded blastocyst stage within 96 hours	Passed	Non-embryotoxic	SPAR-4398-0624
Human Sperm Survival Assay (HSSA)	≥75% sperm motility after 24 hours proven	Passed	Stable human sperm motility	S8029-0624SPAR

Quality Management system according to ISO 13485:2016


GOosafe with SparMED!

Manuela Jagieniak

Quality Control Department

SparMED ApS

Issue date: 18/06/2024



Sparmed ApS

Page 1 of 1

Ryttermarken 2, 3520 Farum, DK +45 39 40 25 03 info@sparmed.dk www.sparmed.dk

Corp. VAT no.: DK 30 89 85 75