

# Certificate of Analysis

**Product name:** Oosafe® Disinfectant for Co2 Incubators and Laminar Flow Hoods

**Product type:** Biocide

**Product code:** OODIH

**Batch/LOT No.:** 60975

**Production date:** 11.2024

**Product expiry date:** 11.2027

TESTING METHOD	ACCEPTANCE CRITERIA	RESULTS	DOCUMENT NO.
Color / Visual	Clear liquid	Clear liquid	QC_4.1-0325
Odor / Olfactive	mild	mild	PRD-F01
ph (2% active substance in water)/ ph meter	7.0-8.5	7.7	PRD-F01
Density, 20 °C, gr/cm <sup>3</sup> / Density measuring device	1.01-1.10	1.07	PRD-F01

TESTING METHOD	ACCEPTANCE CRITERIA	RESULTS	CONCLUSION	DOCUMENT NO.
Mouse Embryo Assay Test (MEA)	Over 80% embryo development to the expanded blastocyst stage within 96 hours	Passed	Non-embryotoxic	SPAR-691-25-03
Human Sperm Survival Assay (HSSA)	≥75% sperm motility after 24 hours proven	Passed	Stable human sperm motility	SPAR-690-25-03

Quality Management system according to ISO 13485:2016

**GOosafe with SparMED!**

Manuela Jagieniak

Quality Control Department

SparMED ApS

**Issue date:** 10.04.2025



**Sparmed ApS**

Page 1 of 1

Ryttermarken 2, 3520 Farum, DK +45 39 40 25 03 info@sparmed.dk www.sparmed.dk

Corp. VAT no.: DK 30 89 85 75