

Certificate of Analysis

Product name: OOSAFE® SURFACE DISINFECTANT

Product code: OODSF

Batch/LOT No.: 70894

Production date: 2019-10

Expiry date: 2022-10

Analysis / Method	Specification / Reference No.	Result	Completion date
Appearance / Visual	Clear, slightly yellow liquid	Clear, slightly yellow liquid	NA
Odor / Olfactive	Slight smell of propionic acid	Slight smell of propionic acid	NA
pH-Rate (2% active substance in water) / pH meter (DIN EN 1262)	5,5-7,0	6,2	NA
Density, 20 °C, gr/cm ³ / Density measuring device	1,001-1,05	1,020	NA
Quaternary Ammonium Compounds / QUAT test stripe		< 7 %	NA
Additives		< 6,0 %	NA
Water		the rest	NA
Proven non-embryotoxic by Mouse Embryo Assay Test. Over 80% embryo development to the expanded blastocyst stage within 96 hours.	SPAR-2887-1119	Passed	15-11-2019
Proved stable human sperm motility: ≥75% sperm motility after 24 hours proven.	S4133-1119SPAR	Passed	13-11-2019

Quality Control according to ISO 13458:2016

GOosafe with SparMED!

Inga Sileikaite,
Quality Control Department
SparMED ApS
Date: 18-Nov-2019

