

Certificate of Analysis

Product name: OOSAFE Surface disinfectant

Product code: OODSF

Batch/LOT No.: 70904

Production date: 2020-05

Expiry date: 2023-05

Analysis / Method	Specification / Reference No.	Result	Completion date
Appearance / Visual	Clear, slightly yellow liquid	Clear, slightly yellow liquid	NA
Odor / Olfactive	Slight smell of propionic acid	Slight smell of propionic acid	NA
pH-Rate (2% active substance in water) / pH meter (DIN EN 1262)	5,5 – 7,0	6,3	NA
Density, 20 °C, gr/cm ³ / Density measuring device	1,001 – 1,05	1,020	NA
Quaternary Ammonium Compounds / QUAT test stripe		< 7,0 %	NA
Additives		< 6,0 %	NA
Water		the rest	NA
Proven non-embryotoxic by Mouse Embryo Assay Test. Over 80% embryo development to the expanded blastocyst stage within 96 hours.	SPAR-4055-0520	Passed	31-05-2020
Proved stable human sperm motility: ≥75% sperm motility after 24 hours proven.	S4471-0520SPAR	Passed	29-05-2020

Quality Control according to ISO 13458:2016

GOosafe with SparMED!

Florencia Garcia
Quality Control Department
SparMED ApS
Date: 02-06-2020


