

Certificate of Analysis

Product name: OOSAFE® LABORATORY SURFACE DISINFECTANT

Product code: OODSF

Batch/LOT No.: 70914

Production date: 2020-11

Expiry date: 2023-11

Analysis / Method	Specification / Reference No.	Result	Completion date
Appearance / Visual	Clear, slightly yellow liquid	Clear, slightly yellow liquid	NA
Odor / Olfactive	Slight smell of propionic acid	Slight smell of propionic acid	NA
pH-Rate (2% active substance in water) / pH meter (DIN EN 1262)	5,5 – 7,0	6,4	NA
Density, 20 °C, gr/cm ³ / Density measuring device	1,01 – 1,15	1,02	NA
Quaternary Ammonium Compounds / QUAT test stripe		< 7,0 %	NA
Additives		< 6,0 %	NA
Water		The rest	NA
Proven non-embryotoxic by Mouse Embryo Assay Test. Over 80% embryo development to the expanded blastocyst stage within 96 hours.	SPAR-5217-0221	Passed	05-02-2021
Proved stable human sperm motility: ≥75% sperm motility after 24 hours proven.	S4837-0221SPAR	Passed	03-02-2021

Quality Control according to ISO 13458:2016

GOosafe with SparMED!

Florencia Garcia
Quality Control Department
SparMED ApS
Date: 08-02-2021


